

# Start the python interpreter and load the module

```
import pynmrstar
```

Load an entry directly from the BMRB database

```
>>> ent15000 = pynmrstar.Entry.from_database(15000)
```

View the heirarchy of saveframes and loops as a tree

```
>>> ent15000.print_tree()
<pynmrstar.Entry '15000' from_database(15000)>
  [0] <pynmrstar.Saveframe 'entry_information'>
    [0] <pynmrstar.Loop '_Entry_author'>
    [1] <pynmrstar.Loop '_SG_project'>
    [2] <pynmrstar.Loop '_Struct_keywords'>
    [3] <pynmrstar.Loop '_Data_set'>
    [4] <pynmrstar.Loop '_Datum'>
    [5] <pynmrstar.Loop '_Release'>
    [6] <pynmrstar.Loop '_Related_entries'>
  [1] <pynmrstar.Saveframe 'citation_1'>
    [0] <pynmrstar.Loop '_Citation_author'>
  [2] <pynmrstar.Saveframe 'assembly'>
    [0] <pynmrstar.Loop '_Entity_assembly'>
  [3] <pynmrstar.Saveframe 'F5-Phe-cVHP'>
```

```
[0] <pymrstar.Loop '_Entity_db_link'>
[1] <pymrstar.Loop '_Entity_comp_index'>
[2] <pymrstar.Loop '_Entity_poly_seq'>
[4] <pymrstar.Saveframe 'natural_source'>
    [0] <pymrstar.Loop '_Entity_natural_src'>
[5] <pymrstar.Saveframe 'experimental_source'>
    [0] <pymrstar.Loop '_Entity_experimental_src'>
[6] <pymrstar.Saveframe 'chem_comp_PHF'>
    [0] <pymrstar.Loop '_Chem_comp_descriptor'>
    [1] <pymrstar.Loop '_Chem_comp_atom'>
    [2] <pymrstar.Loop '_Chem_comp_bond'>
[7] <pymrstar.Saveframe 'unlabeled_sample'>
    [0] <pymrstar.Loop '_Sample_component'>
[8] <pymrstar.Saveframe 'selectively_labeled_sample'>
    [0] <pymrstar.Loop '_Sample_component'>
[9] <pymrstar.Saveframe 'sample_conditions'>
    [0] <pymrstar.Loop '_Sample_condition_variable'>
[10] <pymrstar.Saveframe 'NMRPipe'>
    [0] <pymrstar.Loop '_Vendor'>
    [1] <pymrstar.Loop '_Task'>
[11] <pymrstar.Saveframe 'PIPP'>
    [0] <pymrstar.Loop '_Vendor'>
    [1] <pymrstar.Loop '_Task'>
[12] <pymrstar.Saveframe 'SPARKY'>
    [0] <pymrstar.Loop '_Vendor'>
    [1] <pymrstar.Loop '_Task'>
[13] <pymrstar.Saveframe 'CYANA'>
    [0] <pymrstar.Loop '_Vendor'>
```

```
[1] <pymrstar.Loop '_Task'>
[14] <pymrstar.Saveframe 'X-PLOR_NIH'>
      [0] <pymrstar.Loop '_Vendor'>
            [1] <pymrstar.Loop '_Task'>
                  [15] <pymrstar.Saveframe 'spectrometer_1'>
                        [16] <pymrstar.Saveframe 'spectrometer_2'>
                              [17] <pymrstar.Saveframe 'spectrometer_3'>
                                    [18] <pymrstar.Saveframe 'spectrometer_4'>
                                          [19] <pymrstar.Saveframe 'spectrometer_5'>
                                                [20] <pymrstar.Saveframe 'spectrometer_6'>
                                                      [21] <pymrstar.Saveframe 'NMR_spectrometer_list'>
                                                            [0] <pymrstar.Loop '_NMR_spectrometer_view'>
                                                                  [22] <pymrstar.Saveframe 'experiment_list'>
                                                                        [0] <pymrstar.Loop '_Experiment'>
                                                                              [23] <pymrstar.Saveframe 'chemical_shift_reference_1'>
                                                                                    [0] <pymrstar.Loop '_Chem_shift_ref'>
                                                                                          [24] <pymrstar.Saveframe 'assigned_chem_shift_list_1'>
                                                                                                [0] <pymrstar.Loop '_Chem_shift_experiment'>
                                                                                                      [1] <pymrstar.Loop '_Atom_chem_shift'>
```

There is a shorthand way to access saveframes by name (look at the tree for the saveframe names):

```
>>> ent15000['entry_information']
<pymrstar.Saveframe 'entry_information'>
```

And a shorthand way to access loops by category (only one loop of a given

category can exist in a saveframe):

```
>>> ent15000['entry_information']['_Entry_author']  
<pynmrstar.Loop '_Entry_author'>
```

Next, we'll load the same entry again from the database, remove a saveframe, and compare it to the original.

```
>>> ent15000_2 = pynmrstar.Entry.from_database(15000)  
>>> del ent15000_2['entry_information']  
>>> pynmrstar.diff(ent15000, ent15000_2)  
The number of saveframes in the entries are not equal: '25' v  
s '24'.  
No saveframe with name 'entry_information' in other entry.
```

Let's look at a loop's tags and its data:

```
>>> ent15000['entry_information']['_Entry_author'].get_tag_na  
mes()  
[u'_Entry_author.Ordinal',  
 u'_Entry_author.Given_name',  
 u'_Entry_author.Family_name',  
 u'_Entry_author.First_initial',  
 u'_Entry_author.Middle_initials',  
 u'_Entry_author.Family_title',  
 u'_Entry_author.Entry_ID']
```

Get the first names of the authors using direct saveframe name and loop reference:

```
>>> ent15000['entry_information']['_Entry_author'].get_data_by_tag('Given_name')  
[[u'Claudia', u'Gabriel', u'Erik', u'Samuel', u'John']]
```

Get the first and last names of the authors by providing multiple tags to `get_data_by_tag`.

```
>>> ent15000['entry_information']['_Entry_author'].get_data_by_tag(['Given_name', 'Family_name'])  
[[u'Claudia', u'Cornilescu'],  
 [u'Gabriel', u'Cornilescu'],  
 [u'Erik', u'Hadley'],  
 [u'Samuel', u'Gellman'],  
 [u'John', u'Markley']]
```

Write the modified entry to disk in NMR-STAR format:

```
>>> ent15000_2.write_to_file("example_entry.str")
```

Get a list of validation errors. (The line numbers are only available if an entry is loaded from a file. When an entry is loaded from the API the line numbers are not preserved.)

```
>>> ent15000.validate()  
["Value cannot be NULL but is: '_Chem_comp.Provenance': '.' on  
line 'None'."]
```

Here is how to create a loop from scratch

```
>>> new_loop = pynmrstar.Loop.from_scratch()  
>>> new_loop.add_tag("_Example.age")  
>>> new_loop.add_tag("name")  
>>> new_loop.add_tag("description")  
>>> new_loop.add_tag("notes")
```

Alternatively, you could replace above with:

```
>>> new_loop.add_tag(["_Example.age", "name", "description", "no  
tes"])
```

Now let's add data to the loop

```
>>> new_loop.add_data([29, "Jon", "A BRMB employee", None])
```

Notice that data is automatically encapsulated as necessary to meet the NMR-STAR format (quotes around the data containing a space). You never have to worry about encapsulating data you insert to make it syntactically valid STAR. Notice also that python None types are automatically converted to the NMR-STAR “null” value, “.”.

```
>>> print new_loop
loop_
  _Example.age
  _Example.name
  _Example.description
  _Example.notes

29    Jon    'A BRMB employee'

stop_
```

Add the loop to the entry\_information saveframe

```
>>> ent15000['entry_information'].add_loop(new_loop)
```

See that the loop has been added to the saveframe

```
>>> ent15000['entry_information'].print_tree()
<pymrstar.Saveframe 'entry_information'>
[0] <pymrstar.Loop '_Entry_author'>
[1] <pymrstar.Loop '_SG_project'>
[2] <pymrstar.Loop '_Struct_keywords'>
[3] <pymrstar.Loop '_Data_set'>
[4] <pymrstar.Loop '_Datum'>
[5] <pymrstar.Loop '_Release'>
[6] <pymrstar.Loop '_Related_entries'>
```

```
[7] <pynmrstar.Loop '_Example'>
```

View the value of the tag **ID** in the assembly saveframe

```
>>> ent15000['assembly'].get_tag('ID')  
1
```

To get the NMR-STAR representation of any object, just request its string representation:

```
>>> print ent15000['assembly']  
#####  
# Molecular system (assembly) description #  
#####  
  
save_assembly  
_Assembly.Sf_category          assembly  
_Assembly.Sf_framecode        assembly  
_Assembly.Entry_ID             15000  
_Assembly.ID                   1  
_Assembly.Name                 F5-Phe-cVHP  
_Assembly.BMRB_code            .  
_Assembly.Number_of_components 1  
_Assembly.Organic_ligands      .  
_Assembly.Metal_ions           .  
_Assembly.Non_standard_bonds   .  
_Assembly.Ambiguous_conformational_states .
```



_Assembly.Ambiguous_chem_comp_sites	.
_Assembly.Molecules_in_chemical_exchange	.
_Assembly.Paramagnetic	no
_Assembly.Thiol_state	'all free'
_Assembly.Molecular_mass	.
_Assembly.Enzyme_commission_number	.
_Assembly.Details	.
_Assembly.DB_query_date	.
_Assembly.DB_query_revised_last_date	.

loop\_

_Entity_assembly.ID
_Entity_assembly.Entity_assembly_name
_Entity_assembly.Entity_ID
_Entity_assembly.Entity_label
_Entity_assembly.Asym_ID
_Entity_assembly.PDB_chain_ID
_Entity_assembly.Experimental_data_reported
_Entity_assembly.Physical_state
_Entity_assembly.Conformational_isomer
_Entity_assembly.Chemical_exchange_state
_Entity_assembly.Magnetic_equivalence_group_code
_Entity_assembly.Role
_Entity_assembly.Details
_Entity_assembly.Entry_ID
_Entity_assembly.Assembly_ID

1	F5-Phe-cVHP	1	\$F5-Phe-cVHP	K	.	yes	nativ
---	-------------	---	---------------	---	---	-----	-------

```
e    no    no    .    .    .    15000    1
```

```
stop_
```

```
save_
```

## Reading Spectral Peaks from a NMR-STAR file:

First load the file and get a list of the peak list saveframes

```
>>> ent6577 = pynmrstar.Entry.from_database(6577)
>>> spectral_peaks = ent6577.get_saveframes_by_category('spectral_peak_list')
```

Lets look at how many spectral peak list saveframes we have

```
>>> spectral_peaks
[<pynmrstar.Saveframe 'peak_list_1'>,
 <pynmrstar.Saveframe 'peak_list_2'>,
 <pynmrstar.Saveframe 'peak_list_3'>]
```

For this demo we'll just look at one individual peak list

```
>>> peak1 = spectral_peaks[0]
```

We can see what loops this peak list saveframe contains

```
>>> peak1.print_tree()
<pymrstar.Saveframe 'peak_list_1'>
  [0] <pymrstar.Loop '_Spectral_peak_software'>
  [1] <pymrstar.Loop '_Peak_general_char'>
  [2] <pymrstar.Loop '_Peak_char'>
  [3] <pymrstar.Loop '_Assigned_peak_chem_shift'>
```

Let's see what the `_Peak_char` loop looks like in NMR-STAR format

```
>>> print peak1['_Peak_char']
loop_
  _Peak_char.Peak_ID
  _Peak_char.Spectral_dim_ID
  _Peak_char.Chem_shift_val
  _Peak_char.Chem_shift_val_err
  _Peak_char.Line_width_val
  _Peak_char.Line_width_val_err
  _Peak_char.Phase_val
  _Peak_char.Phase_val_err
  _Peak_char.Decay_rate_val
  _Peak_char.Decay_rate_val_err
  _Peak_char.Coupling_pattern
  _Peak_char.Bounding_box_upper_val
  _Peak_char.Bounding_box_lower_val
  _Peak_char.Bounding_box_range_val
```

```

_Peak_char.Details
_Peak_char.Derivation_method_ID
_Peak_char.Entry_ID
_Peak_char.Spectral_peak_list_ID

```

```

1      1    9.857    .    .    .    .    .    .    .    .    .    .
.      .    .    6577    1
1      2    4.922    .    .    .    .    .    .    .    .    .    .
.      .    .    6577    1
2      1    9.857    .    .    .    .    .    .    .    .    .    .
.      .    .    6577    1
2      2    2.167    .    .    .    .    .    .    .    .    .    .
.      .    .    6577    1
...etc...

```

That is more information than we want right now. Lets get just the columns we need (we'll get a list of lists, each inner list corresponds to a row:

```

>>> our_data = peak1['_Peak_char'].get_data_by_tag(['Peak_ID'
, 'Chem_shift_val'])
>>> print our_data
[[u'1', u'9.857'],
 [u'1', u'4.922'],
 [u'2', u'9.857'],
 [u'2', u'2.167'],
 [u'3', u'9.855'],
 [u'3', u'1.994'],

```

```
...]
```

Excellent! Now we can iterate through each spectral peak and corresponding shift easily. The data is stored as a python list of lists (2 dimensional array) and we can modify or access it any of the normal ways python allows.

```
>>> for x in our_data:
>>>     print "Sprectral chemical shift value is: " + str(x[1]
)
Sprectral chemical shift value is: 9.857
Sprectral chemical shift value is: 4.922
Sprectral chemical shift value is: 9.857
...
```

It is also easy to dump the table in a loop as a CSV

```
>>> print peak1['_Peak_char'].get_data_as_csv()
_Peak_char.Peak_ID,_Peak_char.Spectral_dim_ID,_Peak_char.Chem
_shift_val,_Peak_char.Chem_shift_val_err,_Peak_char.Line_widt
h_val,_Peak_char.Line_width_val_err,_Peak_char.Phase_val,_Pea
k_char.Phase_val_err,_Peak_char.Decay_rate_val,_Peak_char.Dec
ay_rate_val_err,_Peak_char.Coupling_pattern,_Peak_char.Boundi
ng_box_upper_val,_Peak_char.Bounding_box_lower_val,_Peak_char
.Bounding_box_range_val,_Peak_char.Details,_Peak_char.Derivat
ion_method_ID,_Peak_char.Entry_ID,_Peak_char.Spectral_peak_li
st_ID
```

1, 1, 9.857, . . . . ., 6577, 1

1, 2, 4.922, . . . . ., 6577, 1

2, 1, 9.857, . . . . ., 6577, 1

2, 2, 2.167, . . . . ., 6577, 1

3, 1, 9.855, . . . . ., 6577, 1

3, 2, 1.994, . . . . ., 6577, 1

. . .