PyNMRSTAR

A Python module for reading, writing, and manipulating NMR-STAR files.

build passing

Python versions supported: 2.6, 2.7, 3.3, 3.4, 3.5, and 3.6

Overview

This library was developed by the BMRB to give the Python-using NMR community tools to work with the NMR-STAR data format. It is used internally and is actively maintained. The library is thoroughly documented such that calling help(object_or_method) from an interactive python session will print the documentation for the object or method.

That same documentation, as well as some notes on module-level variables is located here.

Finally, there are several command-line based tools developed to enable simple queries to pull data out of an NMR-STAR file. Those tools also serve as great examples of how to use the library. You can view those here.

Introduction to NMR-STAR

To understand how the library works, you first need to understand the NMR-STAR terminology and file format. If you are already familiar with NMR-STAR, feel free to jump ahead to the section on this library.

A NMR-STAR entry/file is componsed of one or more saveframes (conceptually you should think of a saveframe as a data block), each of which contain tags and loops. There can only be one of each tag in a saveframe. If a tag has multiple values, the only way to represent it is to place it inside a loop. A loop is simply a set of tags with multiple values.

Therefore, hierarchically, you can picture a NMR-STAR file as a tree where the entry is the trunk, the large branches are the saveframes, and each saveframe may contain one or more loops - the branches.

Here is a very simple example of a NMR-STAR file:

```
data_dates

save_special_dates_saveframe_1

__Special_Dates.Type Holidays
loop_

__Events.Date

__Events.Desciption

12/31/2017 "New Year's Eve"

01/01/2018 "New Year's Day"

stop_
save_
```

In the previous example, the entry name is dates because that is what follows the data_ tag. Next, there is one saveframe, with a name of special_dates_saveframe_1 and a tag prefix (which corresponds to the saveframe category) of Special Dates. There is one tag in the saveframe,

with a tag name of Type and a value of Holidays. There is also one loop of category events that has information about two different events (though an unlimited number of events could be present).

The first datum in each row corresponds to the first tag, Date, and the second corresponds to the second tag, Description.

Values in NMR-STAR format need to be quoted if they contain a space, tab, vertical tab, or newline in the value. This library takes care of that for you, but it is worth knowing. That is why in the example the dates are not quoted, but the event descriptions are.

Quick Start to PyNMRSTAR

First, pull up an interactive python session and import the module:

```
>>> import pynmrstar
```

There are many ways to load an NMR-STAR entry, but lets focus on the most common two.

From the BMRB API (loads the most up to date version of an entry from the BMRB API):

```
>>> entry15000 = pynmrstar.Entry.from_database(15000)
```

From a file:

```
>>> entry = pynmrstar.Entry.from_file("/location/of/the/file.
str")
```

Continuing on we will assume you have loaded entry 15000 from the API using the from_database command.

Writing out a modified entry or saveframe to file is just as easy:

```
>>> entry15000.write_to_file("output_file_name.str")
```

Viewing the structure of the entry

To see the overall structure of the entry, use the print_tree() method.

```
>>> entry15000.print_tree()
<pynmrstar.Entry '15000' from_database(15000)>
    [0] <pynmrstar.Saveframe 'entry_information'>
        [0] <pynmrstar.Loop '_Entry_author'>
        [1] <pynmrstar.Loop '_SG_project'>
        [2] <pynmrstar.Loop '_Struct_keywords'>
        [3] <pynmrstar.Loop '_Data_set'>
        [4] <pynmrstar.Loop '_Datum'>
        [5] <pynmrstar.Loop '_Release'>
        [6] <pynmrstar.Loop '_Related_entries'>
        [1] <pynmrstar.Saveframe 'citation_1'>
        [0] <pynmrstar.Loop '_Citation_author'>
```

```
[2] <pynmrstar.Saveframe 'assembly'>
    [0] <pynmrstar.Loop ' Entity assembly'>
[3] <pynmrstar.Saveframe 'F5-Phe-cVHP'>
   [0] <pynmrstar.Loop ' Entity db link'>
    [1] <pynmrstar.Loop '_Entity_comp_index'>
    [2] <pynmrstar.Loop ' Entity poly seq'>
[4] <pynmrstar.Saveframe 'natural source'>
    [0] <pynmrstar.Loop '_Entity_natural_src'>
[5] <pynmrstar.Saveframe 'experimental source'>
    [0] <pynmrstar.Loop ' Entity experimental src'>
[6] <pynmrstar.Saveframe 'chem comp PHF'>
   [0] <pynmrstar.Loop ' Chem comp descriptor'>
    [1] <pynmrstar.Loop ' Chem comp atom'>
    [2] <pynmrstar.Loop ' Chem comp bond'>
[7] <pynmrstar.Saveframe 'unlabeled sample'>
   [0] <pynmrstar.Loop ' Sample component'>
[8] <pynmrstar.Saveframe 'selectively labeled sample'>
    [0] <pynmrstar.Loop ' Sample component'>
[9] <pynmrstar.Saveframe 'sample conditions'>
    [0] <pynmrstar.Loop ' Sample condition variable'>
[10] <pynmrstar.Saveframe 'NMRPipe'>
    [0] <pynmrstar.Loop ' Vendor'>
    [1] <pynmrstar.Loop ' Task'>
[11] <pynmrstar.Saveframe 'PIPP'>
    [0] <pynmrstar.Loop ' Vendor'>
    [1] <pynmrstar.Loop ' Task'>
[12] <pynmrstar.Saveframe 'SPARKY'>
    [0] <pynmrstar.Loop ' Vendor'>
```

```
[1] <pynmrstar.Loop ' Task'>
[13] <pynmrstar.Saveframe 'CYANA'>
    [0] <pynmrstar.Loop ' Vendor'>
    [1] <pynmrstar.Loop ' Task'>
[14] <pynmrstar.Saveframe 'X-PLOR NIH'>
    [0] <pynmrstar.Loop ' Vendor'>
    [1] <pynmrstar.Loop ' Task'>
[15] <pynmrstar.Saveframe 'spectrometer 1'>
[16] <pynmrstar.Saveframe 'spectrometer 2'>
[17] <pynmrstar.Saveframe 'spectrometer 3'>
[18] <pynmrstar.Saveframe 'spectrometer 4'>
[19] <pynmrstar.Saveframe 'spectrometer 5'>
[20] <pynmrstar.Saveframe 'spectrometer 6'>
[21] <pynmrstar.Saveframe 'NMR spectrometer list'>
    [0] <pynmrstar.Loop ' NMR spectrometer view'>
[22] <pynmrstar.Saveframe 'experiment list'>
    [0] <pynmrstar.Loop ' Experiment'>
[23] <pynmrstar.Saveframe 'chemical shift reference 1'>
    [0] <pynmrstar.Loop ' Chem shift ref'>
[24] <pynmrstar.Saveframe 'assigned chem shift list 1'>
    [0] <pynmrstar.Loop '_Chem_shift_experiment'>
    [1] <pynmrstar.Loop ' Atom chem shift'>
```

You can see that there are 24 saveframes, and each saveframe contains some number of loops.

Accessing saveframes and loops

There are several ways to access saveframes and loops depending on what you hope to accomplish.

The interactive session way

When playing with the library, debugging, or learning about NMR-STAR you will most likely find the following method most convenient. Note that it is not the correct pattern to use if you want to iterate all of the data in an entry (for reasons that will be explained below).

You can access the saveframes in an entry directly using their *names*. For example, to get a reference to the spectrometer saveframe named spectrometer_1 you can simply do the following:

```
>>> a_spectrometer = entry15000['spectrometer_1']
```

Note that you can see the saveframe names in the tree printout above.

You can do the same for loops within a saveframe, but for loops you must use their tag category (the part before the period) to access them (note that to get to the <code>Vendor</code> loop we first had to go through its parent saveframe, named <code>X-PLOR_NIH</code> (the <code>X-PLOR_NIH</code> saveframe is of the category <code>software</code> - you'll see where you access the category later and why accessing by category is preferrable).

```
>>> explor_nih_vendor = entry15000['X-PLOR_NIH']['_Vendor']
>>> print explor_nih_vendor
```

These shortcuts are there for your convenience when writing code. The reason you shouldn't use them in production code is because the saveframe names - what you use as a reference - can actually have any arbitrary value. They are fairly consistent, and for certain saveframes are always the same, but for other saveframes users can set them to whatever value they want during the deposition. Therefore the much better way to access data is via the *category*. Note that only one saveframe in an entry can have a given name, but multiple saveframes may be of the same category.

The _ prior to the Vendor loop category is to make it clear you want to access the loop and not a saveframe tag with the name Vendor.

The robust (and recommended) way

A better way to access data is via the category of the data you want to read, or by searching for it with a full tag name. Before going into detail, take a look at what one saveframe from the entry above looks like:

```
#####################################
  Computer software used
save X-PLOR NIH
   Software.Sf category
                         software
  Software.Sf framecode X-PLOR NIH
  Software.Entry ID
                         15000
  Software.ID
                         5
   Software.Name
                          'X-PLOR NIH'
  Software.Version
   Software.Details
  loop
     Vendor.Name
      Vendor.Address
     Vendor.Electronic address
     Vendor.Entry ID
     Vendor.Software ID
     'CD Schwieters, JJ Kuszewski, N Tjandra and GM Clore'
       15000
```

```
loop_

_Task.Task
_Task.Entry_ID
_Task.Software_ID

refinement 15000 5
'structure solution' 15000 5

stop_

save_
```

This is a saveframe describing software that was used during an NMR study. You can see from the saveframe tags that the name of this software package is X-PLOR-NIH. You can see from the tag ID that it is the fifth software saveframe in this entry. The category of this saveframe is "software" which you can see in the Sf_category (short for saveframe category) tag.

This saveframe also has two loops, a vendor loop and a task loop. These are loops rather than free tags as a given software package can have more than one vendor and more than one task it performs.

Reading the software packages

The more robust way to access the data in the software saveframes is by

iterating over all of the software saveframes in the entry and pulling out the data we want. To do this for software, we would write the following:

You can see that this method, <code>get_saveframes_by_category</code> returned all of the software saveframes in the entry. Now we can iterate through them to either pull out data, modify data, or remove data. (One note, each loop category - the text before the period in the loop tags - is unique to its parent saveframe. Therefore you will never find a <code>Task</code> loop in a saveframe with a category of anything other than <code>software</code>. Furthermore, a saveframe can only have one loop of a given category. This means that accessing loops within a saveframe using the category notation is robust and will not lead to you missing a loop.)

The following will combine all the task loops in the entry into CSV format.

```
>>> csv_data = ""
>>> for software_sf in software_saveframes:
>>> print_header = True
```

```
# Wrap this in try/catch because it is not gauranteed
>>>
a software saveframe will have a task loop
       try:
>>>
           csv data += software sf[' Task'].get data as csv(h
>>>
eader=print header)
           print header = False
>>>
       except KeyError:
>>>
           continue
>>>
>>> csv data
' Task.Task, Task.Entry ID, Task.Software ID\nprocessing,1500
0,1\n Task.Task, Task.Entry ID, Task.Software ID\nchemical sh
ift assignment, 15000, 2\ndata analysis, 15000, 2\npeak picking, 1
5000,2\n Task.Task, Task.Entry ID, Task.Software ID\nchemical
 shift assignment, 15000, 3\n Task. Task, Task. Entry ID, Task. So
ftware ID\nstructure solution, 15000, 4\n Task. Task, Task. Entry
ID, Task.Software ID\nrefinement,15000,5\nstructure solution
,15000,5\n'
```

Using get_tag to pull tags directly from an entry

Another way to access data in by using the full tag name. Keep in mind that a full tag contains a category first, a period, and then a tag name. So if we wanted to see all of the various _Task.Task that the software packages associated with this entry performed, a simple way to do so is with the get_tag() method of the entry:

```
>>> entry15000.get_tag('Task.Task')
```

```
[u'processing',
   u'chemical shift assignment',
   u'data analysis',
   u'peak picking',
   u'chemical shift assignment',
   u'structure solution',
   u'refinement',
   u'structure solution']
```

Or to get all of the spectrometer information - get_tags() accepts a list of tags to fetch and returns a dictionary pointing to all the values of each tag, with the order preserved:

```
'_NMR_spectrometer.Model': [u'Avance',
    u'Avance',
    u'Avance',
    u'INOVA',
    u'INOVA',
    u'INOVA']}
```

To view all of the tags in the NMR-STAR schema and their meanings, please go here.

Assigned Chemical Shifts

"I just want to get the chemical shift data as an array - how do I do that?"

Keep in mind that an entry may have multiple sets of assigned chemical shifts. (For examples, there made be two sets of assignments that were made under two different sample conditions.) So to get the chemical shifts it is best to iterate through all the assigned chemical shift loops:

```
>>> cs_result_sets = []
>>> for chemical_shift_loop in entry15000.get_loops_by_catego
ry("Atom_chem_shift"):
>>> cs_result_sets.append(chemical_shift_loop.get_tag(['C
omp_index_ID', 'Comp_ID', 'Atom_ID', 'Atom_type', 'Val', 'Val
_err']))
>>> cs_result_sets
[[[u'2', u'SER', u'H', u'H', u'9.3070', u'0.01'],
```

```
[u'2', u'SER', u'HA', u'H', u'4.5970', u'0.01'],
[u'2', u'SER', u'HB2', u'H', u'4.3010', u'0.01'],
[u'2', u'SER', u'HB3', u'H', u'4.0550', u'0.01'],
[u'2', u'SER', u'CB', u'C', u'64.6000', u'0.1'],
[u'2', u'SER', u'N', u'N', u'121.5800', u'0.1'],
[u'3', u'ASP', u'H', u'H', u'8.0740', u'0.01'],
[u'3', u'ASP', u'HA', u'H', u'4.5580', u'0.01'],
[u'3', u'ASP', u'HB2', u'H', u'2.835', u'0.01'],
...
```

Note that we used the get_tag() method of the loop to only pull out the tags we were concerned with. get_tag() accepts an array of tags in addition to a single tag. The full assigned chemical saveframe loop will contain extra tags you may not need. For example:

```
Atom chem shift. Atom isotope number
      Atom chem shift.Val
      Atom chem shift.Val err
      Atom chem shift. Assign fig of merit
      Atom chem shift. Ambiguity code
      Atom chem shift.Occupancy
      Atom chem shift.Resonance ID
      Atom chem shift. Auth entity assembly ID
      Atom chem shift. Auth asym ID
      Atom chem shift. Auth seg ID
      Atom chem shift. Auth comp ID
      Atom chem shift. Auth atom ID
      Atom chem shift.Details
      Atom chem shift. Entry ID
      _Atom_chem_shift.Assigned chem shift list ID
               1
                   1
                             2
                                  SER
                                        Н
                                               Н
                                                   1
                                                         9.3070
     1
     0.01
                                 . 2
                                          SER
                                               Н
                                                           1500
0
    1
                                               Н
                                                         4.5970
                                  SER
                                        HA
                                                   1
     2
                             2
                                     2
     0.01
                                          SER
                                                HA
                                                           1500
0
    1
                                                         4.3010
                   1
                                  SER
     3
               1
                       2
                             2
                                        HB2
                                               Н
                                                   1
     0.01
                                               HB2
                                                           1500
                                     2
                                          SER
0
    1
```

That is easy to do. When you first load an entry it is by default loaded with all values as strings. To instead load it such that the values match the schema, simply turn on CONVERT DATATYPES prior to loading it.

```
>>> pynmrstar.CONVERT DATATYPES = True
>>> ent15000 = pynmrstar.Entry.from database(15000)
>>> cs result sets = []
>>> for chemical shift loop in entry15000.get loops by catego
ry("Atom chem shift"):
>>> cs result sets.append(chemical shift loop.get tag(['C
omp index ID', 'Comp ID', 'Atom ID', 'Atom type', 'Val', 'Val
err']))
>>> cs result sets
[[[2, u'SER', u'H', u'H', Decimal('9.3070'), Decimal('0.01')]
 [2, u'SER', u'HA', u'H', Decimal('4.5970'), Decimal('0.01')
],
 [2, u'SER', u'HB2', u'H', Decimal('4.3010'), Decimal('0.01'
)],
 [2, u'SER', u'HB3', u'H', Decimal('4.0550'), Decimal('0.01'
)],
  [2, u'SER', u'CB', u'C', Decimal('64.6000'), Decimal('0.1')
],
  [2, u'SER', u'N', u'N', Decimal('121.5800'), Decimal('0.1')
],
```

```
[3, u'ASP', u'H', u'H', Decimal('8.0740'), Decimal('0.01')]

,

[3, u'ASP', u'HA', u'H', Decimal('4.5580'), Decimal('0.01')]
],

[3, u'ASP', u'HB2', u'H', Decimal('2.835'), Decimal('0.01')]
],

[3, u'ASP', u'HB3', u'H', Decimal('2.754'), Decimal('0.01')]
],

[3, u'ASP', u'CA', u'C', Decimal('57.6400'), Decimal('0.1')]
],

[3, u'ASP', u'N', u'N', Decimal('121.1040'), Decimal('0.1')]
],

...
```

This is a great opportunity to point out that if all you want is the chemical shifts, or one or two tags, you may find it significantly easier to use the BMRB API (chemical shift endpoint) to fetch that data directly and ondemand rather than dealing directly with NMR-STAR at all.

Creating new loops and saveframes

This tutorial has so far focused on how to read and access data. This section will focus on how to create new loop and saveframe objects.

Loops

```
There are five ways to make a new loop: <code>from_file()</code>, <code>from_json()</code>, <code>from_scratch()</code>, <code>from_string()</code>, and <code>from_template()</code>. All of these are classmethods. <code>from_scratch()</code> makes a new loop, <code>from_string()</code> parses an NMR-STAR loop from a python string containing NMR-STAR data, <code>from_json()</code> parses a JSON object (reversely, <code>get_json()</code> will get a JSON representation of the loop), <code>from_scratch()</code> makes a completely empty loop, and <code>from_template()</code> makes a loop with the tags prefilled from the BMRB schema based on the provided category. <code>from_file</code>, <code>from_json</code>, and <code>from_string</code> are fairly self-explanatory - see the full documentation if needed for usage.
```

from_scratch()

```
>>> lp = pynmrstar.Loop.from scratch()
>>> print lp
   loop
   stop
>>> lp.add tag(['loop category.tag1', 'loop category.tag2',
loop category.tag3'])
>>> print lp
   loop
      loop category.tag1
      loop category.tag2
      loop category.tag3
```

```
stop_
# Note that when calling add data the length of the array mus
t match the number of tags in the loop
>>> lp.add_data(['value_1', 2, 'value 3'])
>>> print lp
   loop_
     loop category.tag1
     loop category.tag2
      loop category.tag3
     value_1 2 'value 3'
   stop
# Alternatively, you can (with caution) directly modify the a
rray corresponding to the loop data
>>> lp.data = [[1,2,3],[4,5,6]]
>>> print lp
  loop
      loop category.tag1
     loop category.tag2
     _loop_category.tag3
     4 5 6
```

```
stop
```

Note that the loop category was set automatically when the tag

loop_category.tag1 was added. You could have also provided the tag

when creating the loop by providing it as an argument to the optional

category argument to the constructor.

from_template()

This method will create a new loop ready for data with the tags from the BMRB schema corresponding to that loop category.

```
>>> chemical shifts = pynmrstar.Loop.from template('atom chem
 shift list')
>>> print chemical shifts
   loop
      Atom chem shift.ID
      Atom chem shift. Assembly atom ID
      Atom chem shift. Entity assembly ID
      Atom chem shift. Entity ID
      Atom chem shift.Comp index ID
      Atom chem shift.Seq ID
      Atom chem shift.Comp ID
      Atom chem shift.Atom ID
      Atom chem shift. Atom type
      Atom chem shift. Atom isotope number
      Atom chem shift.Val
```

```
Atom chem shift.Val err
  Atom chem shift. Assign fig of merit
  Atom chem shift. Ambiguity code
  Atom chem shift. Ambiguity set ID
  Atom chem shift.Occupancy
  Atom chem shift.Resonance ID
  Atom chem shift. Auth entity assembly ID
   Atom chem shift. Auth asym ID
  Atom chem shift. Auth seq ID
  Atom chem shift. Auth comp ID
  Atom chem shift. Auth atom ID
   Atom chem shift.Details
  Atom chem shift. Entry ID
  Atom chem shift. Assigned chem shift list ID
stop
```

Saveframes

```
There are five ways to make a new loop: from_file(), from_json(),
from_scratch(), from_string(), and from_template(). All of these
are classmethods. from_scratch() makes a new saveframe,
from_string() parses an NMR-STAR saveframe from a python string
containing NMR-STAR data, from_json() parses a JSON object (reversely,
get_json() will get a JSON representation of the saveframe),
from_scratch() makes a completely empty saveframe, and
```

from_template() makes a saveframe with the tags prefilled from the
BMRB schema based on the provided category. from_file, from_json,
and from_string are fairly self-explanatory - see the full documentation if
needed for usage.

from_scratch()

```
# You must provide the saveframe name (the value that comes a
fter "save " at the start of the saveframe and saveframe tag
prefix (the value before the "." in a tag name) when creating
 a saveframe this way
>>> my sf = pynmrstar.Saveframe.from scratch("sf name", "exam
ple sf category")
>>> print my sf
save sf name
save
# Add a tag using the add tag() method. Update=True will over
ride existing tag with the same name. Update=False will raise
 an exception if the tag already exists
>>> my sf.add tag("tagName1", "tagValue1")
>>> print my sf
save sf name
   example sf category.tagName1 tagValue1
save
```

```
>>> my sf.add tag("tagName1", "tagValue2", update=False)
ValueError: There is already a tag with the name 'tagName1'.
>>> my_sf.add_tag("tagName1", "tagValue2", update=True)
>>> print my sf
save sf name
  example sf category.tagName1 tagValue1
save
# Alternatively, you can access or write tag values using dir
ect subset access:
>>> my sf['tagName1']
['tagValue2']
>>> my_sf['tagName2'] = "some value"
>>> print my sf
save sf name
 example sf category.tagName1 tagValue2
   example sf category.tagName2 'some value'
save
# Now add the loop we created before
>>> my sf.add loop(lp)
>>> print my sf
save sf name
  example sf category.tagName1 tagValue2
   _example_sf_category.tagName2 'some value'
   loop
```

```
loop category.tag1
      loop category.tag2
      loop category.tag3
     1 2 3
     4 5 6
   stop
save
# Now write out our saveframe to a file. Optionally specify f
ormat="json" to write in JSON format.
>>> my sf.write to file("file name.str")
>>> my sf.write to file("file name.json", format ="json")
```

from_template()

```
######
        Chemical Shift Ambiguity Index Value Definitions
     #
#
# The values other than 1 are used for those atoms with diffe
rent #
# chemical shifts that cannot be assigned to stereospecific a
toms #
# or to specific residues or chains.
     #
#
    Index Value
                            Definition
#
     #
                     Unique (including isolated methyl proton
#
       1
S,
                           geminal atoms, and geminal methyl
#
                           groups with identical chemical shif
#
ts) #
                           (e.g. ILE HD11, HD12, HD13 protons)
#
     #
       2
                     Ambiguity of geminal atoms or geminal me
#
thyl #
                           proton groups (e.g. ASP HB2 and HB3
```

```
protons, LEU CD1 and CD2 carbons, o
#
     #
                           LEU HD11, HD12, HD13 and HD21, HD22
#
                          HD23 methyl protons)
#
                     Aromatic atoms on opposite sides of
       3
#
                           symmetrical rings (e.g. TYR HE1 and
#
HE2 #
                           protons)
#
     #
                     Intraresidue ambiguities (e.g. LYS HG an
       4
d
    #
                          HD protons or TRP HZ2 and HZ3 proto
ns) #
                     Interresidue ambiguities (LYS 12 vs. LYS
27) #
                     Intermolecular ambiguities (e.g. ASP 31
       6
CA #
                           in monomer 1 and ASP 31 CA in monom
#
er 2 #
                           of an asymmetrical homodimer, duple
X
                           DNA assignments, or other assignmen
#
ts #
                           that may apply to atoms in one or m
```

```
ore #
                        molecule in the molecular assembly)
#
    #
      9
                   Ambiguous, specific ambiguity not define
#
d
#
    #
######
save assigned chemical shifts
   Assigned chem shift list.Sf category
                                                       ass
igned chemical shifts
   _Assigned_chem_shift list.Sf framecode
                                                       ass
igned chemical shifts
   Assigned chem shift list. Entry ID
   Assigned chem shift list.ID
   Assigned chem shift list. Sample condition list ID
   Assigned chem shift list. Sample condition list label
   _Assigned_chem_shift list.Chem shift reference ID
   Assigned chem shift list.Chem shift reference label
   Assigned chem shift list.Chem shift 1H err
   Assigned chem shift list.Chem shift 13C err
   Assigned chem shift list.Chem shift 15N err
   Assigned chem shift list.Chem shift 31P err
   Assigned chem shift list.Chem shift 2H err
   Assigned chem shift list.Chem shift 19F err
   Assigned chem shift list. Error derivation method
```

```
Assigned chem shift list.Details
   Assigned chem shift list.Text data format
  Assigned chem shift list. Text data
  loop
      _Chem_shift_experiment.Experiment ID
      Chem shift experiment. Experiment name
      Chem shift experiment. Sample ID
      Chem shift experiment. Sample label
      Chem shift experiment. Sample state
      Chem shift experiment. Entry ID
      Chem shift experiment. Assigned chem shift list ID
  stop
  loop
      Systematic chem shift offset. Type
      Systematic chem shift offset. Atom type
      Systematic chem shift offset. Atom isotope number
      Systematic chem shift offset.Val
      Systematic chem shift offset.Val err
      Systematic chem shift offset.Entry ID
      Systematic chem shift offset. Assigned chem shift list
ID
  stop
```

```
loop_
   Chem shift software.Software ID
   Chem shift software.Software label
   Chem shift software.Method ID
   Chem shift software. Method label
   Chem shift software. Entry ID
   Chem shift software. Assigned chem shift list ID
stop
loop
   Atom chem shift.ID
   Atom chem shift. Assembly atom ID
   Atom chem shift. Entity assembly ID
   Atom chem shift. Entity ID
   Atom chem shift.Comp index ID
   Atom chem shift. Seq ID
   Atom chem shift.Comp ID
   Atom chem shift.Atom ID
   Atom chem shift. Atom type
   Atom chem shift. Atom isotope number
   Atom chem shift.Val
   Atom chem shift.Val err
   Atom chem shift. Assign fig of merit
   Atom chem shift. Ambiguity code
   Atom chem shift. Ambiguity set ID
```

```
Atom chem shift.Occupancy
      Atom chem shift.Resonance ID
      Atom chem shift. Auth entity assembly ID
      Atom chem shift. Auth asym ID
      Atom chem shift. Auth seq ID
      Atom chem shift. Auth comp ID
      Atom chem shift. Auth atom ID
      Atom chem shift.Details
      Atom chem shift. Entry ID
      Atom chem shift. Assigned chem shift list ID
   stop
   loop
      Ambiguous atom chem shift. Ambiguous shift set ID
      Ambiguous atom chem shift. Atom chem shift ID
      Ambiguous atom chem shift. Entry ID
      _Ambiguous_atom_chem_shift.Assigned chem shift list ID
   stop
save
```

Schema methods

The library makes it easy to add missing tags, sort the tags according to the BMRB schema, and validate the data against the schema. Let's do a simple example of creating a chemical shift loop, adding any missing tags, ordering the tags in the standard order (not required), and then checking for errors.

```
# Create the loop with the proper category
>>> my cs loop = pynmrstar.Loop.from scratch("Atom chem shift
")
# Add the tags we will fill
>>> my cs loop.add tag(['Comp ID', 'Atom ID', 'Comp index ID'
, 'Atom type', 'Val', 'Val err'])
   loop
      Atom chem shift.Comp ID
      Atom chem shift.Atom ID
      Atom chem shift.Comp Index ID
      Atom chem shift. Atom type
      Atom chem shift.Val
      Atom chem shift. Val err
   stop
# Populate the data array
>>> my cs loop.data = [['SER', 'H', '2', 'H', '9.3070', '0.0
1'],
                       ['SER', 'HA', '2', 'H', '4.5970', '0.0
1'],
```

```
['SER', 'HB2', '2', 'H', '4.3010', '0.
01']]
>>> print my cs loop
   loop
     Atom chem shift.Comp ID
      Atom chem shift.Atom ID
     Atom chem shift.Comp Index ID
     Atom chem shift. Atom type
     Atom chem shift.Val
     Atom chem shift. Val err
                2
                    H 9.3070 0.01
     SER
          Н
     SER
          HA
                2
                    H 4.5970 0.01
                2 H
     SER
                      4.3010 0.01
         HB2
   stop
# Now lets sort the tags to match the BMRB schema
>>> my cs loop.sort tags()
# You can see that the Comp index ID tag has been moved to th
e front to match the BMRB standard
>>> print my cs loop
   loop
     Atom chem shift.Comp index ID
     Atom chem shift.Comp ID
     Atom chem shift.Atom ID
     Atom chem shift. Atom type
     Atom chem shift.Val
```

```
Atom chem shift. Val err
        SER
              Н
                    Н
                        9.3070
                                0.01
     2
     2
        SER
             HA
                    Н
                        4.5970 0.01
        SER HB2
                    Н
                        4.3010 0.01
   stop
# Check for any errors - returns a list of errors. No errors
here:
>>> print my cs loop.validate()
# Let us now set 'Comp index ID' to have an invalid value
>>> my cs loop.data[0][0] = "invalid"
# You can see that there is now a validation error - the data
 doesn't match the specified type
>>> print my cs loop.validate()
["Value does not match specification: ' Atom chem shift.Comp
index ID':'invalid' on line '0 tag 0 of loop'.\n
Type spe
cified: int\n
              Regular expression for type: '-?[0-9]+'"]
# If you use the pynmrstar.validate(object) function, it will
 print the report in a human-readable format
>>> pynmrstar.validate(my cs loop)
1: Value does not match specification: ' Atom chem shift.Comp
index ID': 'invalid' on line '0 tag 0 of loop'.
    Type specified: int
     Regular expression for type: '-?[0-9]+'
```

```
# Finally, add in any tags that you didn't have a value for
>>> my cs loop.add missing tags()
# You can see that all the standard "Atom chem shift" loop ta
gs have been added, and their values all set to a logical nul
l value - "."
>>> print my cs loop
   loop
      Atom chem shift.ID
      Atom chem shift. Assembly atom ID
      Atom chem shift. Entity assembly ID
      Atom chem shift. Entity ID
      Atom chem shift.Comp index ID
      Atom chem shift.Seq ID
      Atom chem shift.Comp ID
      Atom chem shift.Atom ID
      Atom chem shift. Atom type
      Atom chem shift.Atom isotope number
      Atom chem shift.Val
      Atom chem shift. Val err
      Atom chem shift. Assign fig of merit
      Atom chem shift. Ambiguity code
      Atom chem shift. Ambiguity set ID
      Atom chem shift.Occupancy
      Atom chem shift.Resonance ID
      Atom chem shift. Auth entity assembly ID
      Atom chem shift. Auth asym ID
      Atom chem shift. Auth seq ID
```

```
Atom chem shift. Auth comp ID
   Atom chem shift. Auth atom ID
   Atom chem shift.Details
   Atom chem shift. Entry ID
   Atom chem shift. Assigned chem shift list ID
                  invalid
                                 SER
                                                      9.3070
                                       Н
                                             Н
0.01
                                 SER
                                       HA
                                                      4.5970
                  2
                                             Н
0.01
                                 SER
                                       HB2
                                                      4.3010
0.01
stop
```

For more examples of PyNMRSTAR library usage, please look here.

For the full documentation of all available methods and classes, please look here.

For any questions or suggestions, please create an issue on the GitHub page.