

PyNMRSTAR

A Python module for reading, writing, and manipulating NMR-STAR files.

build passing

Python versions supported: 2.6, 2.7, 3.3, 3.4, 3.5, and 3.6

Overview

This library was developed by the BMRB to give the Python-using NMR community tools to work with the NMR-STAR data format. It is used internally and is actively maintained. The library is thoroughly documented such that calling `help(object_or_method)` from an interactive python session will print the documentation for the object or method.

That same documentation, as well as some notes on module-level variables is located [here](#).

Finally, there are several command-line based tools developed to enable simple queries to pull data out of an NMR-STAR file. Those tools also serve as great examples of how to use the library. You can view those [here](#).

Introduction to NMR-STAR

To understand how the library works, you first need to understand the NMR-STAR terminology and file format. If you are already familiar with NMR-STAR, feel free to [jump ahead](#) to the section on this library.

A NMR-STAR entry/file is composed of one or more saveframes (conceptually you should think of a saveframe as a data block), each of which contain tags and loops. There can only be one of each tag in a saveframe. If a tag has multiple values, the only way to represent it is to place it inside a loop. A loop is simply a set of tags with multiple values.

Therefore, hierarchically, you can picture a NMR-STAR file as a tree where the entry is the trunk, the large branches are the saveframes, and each saveframe may contain one or more loops - the branches.

Here is a very simple example of a NMR-STAR file:

data_dates	
save_special_dates_saveframe_1	
_Special_Dates.Type	Holidays
loop_	
_Events.Date	
_Events.Description	
12/31/2017 "New Year's Eve"	
01/01/2018 "New Year's Day"	
stop_	
save_	

In the previous example, the entry name is `dates` because that is what follows the `data_` tag. Next, there is one saveframe, with a name of `special_dates_saveframe_1` and a tag prefix (which corresponds to the saveframe category) of `Special_Dates`. There is one tag in the saveframe,

with a tag name of `Type` and a value of `Holidays`. There is also one loop of category `events` that has information about two different events (though an unlimited number of events could be present).

The first datum in each row corresponds to the first tag, `Date`, and the second corresponds to the second tag, `Description`.

Values in NMR-STAR format need to be quoted if they contain a space, tab, vertical tab, or newline in the value. This library takes care of that for you, but it is worth knowing. That is why in the example the dates are not quoted, but the event descriptions are.

Quick Start to PyNMRSTAR

First, pull up an interactive python session and import the module:

```
>>> import pynmrstar
```

There are many ways to load an NMR-STAR entry, but lets focus on the most common two.

From the BMRB API (loads the most up to date version of an entry from the BMRB API):

```
>>> entry15000 = pynmrstar.Entry.from_database(15000)
```

From a file:

```
>>> entry = pynmrstar.Entry.from_file("/location/of/the/file.  
str")
```

Continuing on we will assume you have loaded entry 15000 from the API using the `from_database` command.

Viewing the structure of the entry

To see the overall structure of the entry, use the `print_tree()` method.

```
>>> entry15000.print_tree()  
<pynmrstar.Entry '15000' from_database(15000)>  
  [0] <pynmrstar.Saveframe 'entry_information'>  
    [0] <pynmrstar.Loop '_Entry_author'>  
    [1] <pynmrstar.Loop '_SG_project'>  
    [2] <pynmrstar.Loop '_Struct_keywords'>  
    [3] <pynmrstar.Loop '_Data_set'>  
    [4] <pynmrstar.Loop '_Datum'>  
    [5] <pynmrstar.Loop '_Release'>  
    [6] <pynmrstar.Loop '_Related_entries'>  
  [1] <pynmrstar.Saveframe 'citation_1'>  
    [0] <pynmrstar.Loop '_Citation_author'>  
  [2] <pynmrstar.Saveframe 'assembly'>  
    [0] <pynmrstar.Loop '_Entity_assembly'>  
  [3] <pynmrstar.Saveframe 'F5-Phe-cVHP'>  
    [0] <pynmrstar.Loop '_Entity_db_link'>  
    [1] <pynmrstar.Loop '_Entity_comp_index'>
```

```
[2] <pynmrstar.Loop '_Entity_poly_seq'>
[4] <pynmrstar.Saveframe 'natural_source'>
    [0] <pynmrstar.Loop '_Entity_natural_src'>
[5] <pynmrstar.Saveframe 'experimental_source'>
    [0] <pynmrstar.Loop '_Entity_experimental_src'>
[6] <pynmrstar.Saveframe 'chem_comp_PHF'>
    [0] <pynmrstar.Loop '_Chem_comp_descriptor'>
    [1] <pynmrstar.Loop '_Chem_comp_atom'>
    [2] <pynmrstar.Loop '_Chem_comp_bond'>
[7] <pynmrstar.Saveframe 'unlabeled_sample'>
    [0] <pynmrstar.Loop '_Sample_component'>
[8] <pynmrstar.Saveframe 'selectively_labeled_sample'>
    [0] <pynmrstar.Loop '_Sample_component'>
[9] <pynmrstar.Saveframe 'sample_conditions'>
    [0] <pynmrstar.Loop '_Sample_condition_variable'>
[10] <pynmrstar.Saveframe 'NMRPipe'>
    [0] <pynmrstar.Loop '_Vendor'>
    [1] <pynmrstar.Loop '_Task'>
[11] <pynmrstar.Saveframe 'PIPP'>
    [0] <pynmrstar.Loop '_Vendor'>
    [1] <pynmrstar.Loop '_Task'>
[12] <pynmrstar.Saveframe 'SPARKY'>
    [0] <pynmrstar.Loop '_Vendor'>
    [1] <pynmrstar.Loop '_Task'>
[13] <pynmrstar.Saveframe 'CYANA'>
    [0] <pynmrstar.Loop '_Vendor'>
    [1] <pynmrstar.Loop '_Task'>
[14] <pynmrstar.Saveframe 'X-PLOR_NIH'>
```

```
[0] <pynmrstar.Loop '_Vendor'>
[1] <pynmrstar.Loop '_Task'>
[15] <pynmrstar.Saveframe 'spectrometer_1'>
[16] <pynmrstar.Saveframe 'spectrometer_2'>
[17] <pynmrstar.Saveframe 'spectrometer_3'>
[18] <pynmrstar.Saveframe 'spectrometer_4'>
[19] <pynmrstar.Saveframe 'spectrometer_5'>
[20] <pynmrstar.Saveframe 'spectrometer_6'>
[21] <pynmrstar.Saveframe 'NMR_spectrometer_list'>
      [0] <pynmrstar.Loop '_NMR_spectrometer_view'>
[22] <pynmrstar.Saveframe 'experiment_list'>
      [0] <pynmrstar.Loop '_Experiment'>
[23] <pynmrstar.Saveframe 'chemical_shift_reference_1'>
      [0] <pynmrstar.Loop '_Chem_shift_ref'>
[24] <pynmrstar.Saveframe 'assigned_chem_shift_list_1'>
      [0] <pynmrstar.Loop '_Chem_shift_experiment'>
      [1] <pynmrstar.Loop '_Atom_chem_shift'>
```

You can see that there are 24 saveframes, and each saveframe contains some number of loops.

Accessing saveframes and loops

There are several ways to access saveframes and loops depending on what you hope to accomplish.

The interactive session way

When playing with the library, debugging, or learning about NMR-STAR you will most likely find the following method most convenient. Note that it is not the correct pattern to use if you want to iterate all of the data in an entry (for reasons that will be explained below).

You can access the saveframes in an entry directly using their *names*. For example, to get a reference to the spectrometer saveframe named `spectrometer_1` you can simply do the following:

```
>>> a_spectrometer = entry15000['spectrometer_1']
```

Note that you can see the saveframe names in the tree printout above.

You can do the same for loops within a saveframe, but for loops you must use their tag category (the part before the period) to access them (note that to get to the `Vendor` loop we first had to go through its parent saveframe, named `X-PLOR_NIH` (the `X-PLOR_NIH` saveframe is of the category `software` - you'll see where you access the category later and why accessing by category is preferable).

```
>>> explor_nih_vendor = entry15000['X-PLOR_NIH']['_Vendor']
>>> print explor_nih_vendor
loop_
  _Vendor.Name
  _Vendor.Address
  _Vendor.Electronic_address
  _Vendor.Entry_ID
```

```
_Vendor.Software_ID
```

```
'CD Schwieters, JJ Kuszewski, N Tjandra and GM Clore'
```

```
. . 15000 5
```

```
stop_
```

These shortcuts are there for your convenience when writing code. The reason you shouldn't use them in production code is because the saveframe names - what you use as a reference - can actually have any arbitrary value. They are fairly consistent, and for certain saveframes are always the same, but for other saveframes users can set them to whatever value they want during the deposition. Therefore the much better way to access data is via the *category*. Note that only one saveframe in an entry can have a given name, but multiple saveframes may be of the same category.

The `_` prior to the `Vendor` loop category is to make it clear you want to access the loop and not a saveframe tag with the name `Vendor`.

The robust (and recommended) way

A better way to access data is via the category of the data you want to read, or by searching for it with a full tag name. Before going into detail, take a look at what one saveframe from the entry above looks like:

```
#####
```



```
# Computer software used #
```

```
#####
```

```
save_X-PLOR_NIH
```

```
_Software.Sf_category    software
```

```
_Software.Sf_framecode   X-PLOR_NIH
```

```
_Software.Entry_ID       15000
```

```
_Software.ID             5
```

```
_Software.Name           'X-PLOR NIH'
```

```
_Software.Version        .
```

```
_Software.Details        .
```

```
loop_
```

```
_Vendor.Name
```

```
_Vendor.Address
```

```
_Vendor.Electronic_address
```

```
_Vendor.Entry_ID
```

```
_Vendor.Software_ID
```

```
'CD Schwieters, JJ Kuszewski, N Tjandra and GM Clore'
```

```
. . 15000 5
```

```
stop_
```

```
loop_
```

```
_Task.Task
```

```
_Task.Entry_ID
```

```
_Task.Software_ID
```

refinement	15000	5
'structure solution'	15000	5
stop_		
save_		

This is a saveframe describing software that was used during an NMR study. You can see from the saveframe tags that the name of this software package is X-PLOR-NIH. You can see from the tag **ID** that it is the fifth software saveframe in this entry. The category of this saveframe is “software” which you can see in the **Sf_category** (short for saveframe category) tag.

This saveframe also has two loops, a vendor loop and a task loop. These are loops rather than free tags as a given software package can have more than one vendor and more than one task it performs.

Reading the software packages

The more robust way to access the data in the software saveframes is by iterating over all of the software saveframes in the entry and pulling out the data we want. To do this for software, we would write the following:

```
software_saveframes = entry15000.get_saveframes_by_category('
software')
software_saveframes
```

```
[<pynmrstar.Saveframe 'NMRPipe'>,  
 <pynmrstar.Saveframe 'PIPP'>,  
 <pynmrstar.Saveframe 'SPARKY'>,  
 <pynmrstar.Saveframe 'CYANA'>,  
 <pynmrstar.Saveframe 'X-PLOR_NIH'>]
```

You can see that this method, `get_saveframes_by_category` returned all of the software saveframes in the entry. Now we can iterate through them to either pull out data, modify data, or remove data. (One note, each loop category - the text before the period in the loop tags - is unique to its parent saveframe. Therefore you will never find a `Task` loop in a saveframe with a category of anything other than `software`. Furthermore, a saveframe can only have one loop of a given category. This means that accessing loops within a saveframe using the category notation is robust and will not lead to you missing a loop.)

The following will combine all the task loops in the entry into CSV format.

```
csv_data = ""  
for software_sf in software_saveframes:  
    print_header = True  
    # Wrap this in try/catch because it is not guaranteed a s  
software saveframe will have a task loop  
    try:  
        csv_data += software_sf['_Task'].get_data_as_csv(head  
er=print_header)  
        print_header = False
```

```
except KeyError:
```

```
    continue
```

```
csv_data
```

```
'_Task.Task,_Task.Entry_ID,_Task.Software_ID\nprocessing,15000,1\n_Task.Task,_Task.Entry_ID,_Task.Software_ID\nchemical shift assignment,15000,2\ndata analysis,15000,2\npeak picking,15000,2\n_Task.Task,_Task.Entry_ID,_Task.Software_ID\nchemical shift assignment,15000,3\n_Task.Task,_Task.Entry_ID,_Task.Software_ID\nstructure solution,15000,4\n_Task.Task,_Task.Entry_ID,_Task.Software_ID\nrefinement,15000,5\nstructure solution,15000,5\n'
```

Using get_tag to pull tags directly from an entry

Another way to access data in by using the full tag name. Keep in mind that a full tag contains a category first, a period, and then a tag name. So if we wanted to see all of the various `_Task.Task` that the software packages associated with this entry performed, a simple way to do so is with the `get_tag()` method of the entry:

```
print entry15000.get_tag('Task.Task')
```

```
[u'processing',  
 u'chemical shift assignment',  
 u'data analysis',  
 u'peak picking',  
 u'chemical shift assignment',
```

```
u'structure solution',  
u'refinement',  
u'structure solution']
```

Or to get all of the spectrometer information - `get_tags()` accepts a list of tags to fetch and returns a dictionary pointing to all the values of each tag, with the order preserved:

```
entry15000.get_tags(['_NMR_spectrometer.Manufacturer', '_NMR_  
spectrometer.Model', '_NMR_spectrometer.Field_strength'])  
{'_NMR_spectrometer.Field_strength': [u'500',  
u'500',  
u'750',  
u'600',  
u'800',  
u'900'],  
'_NMR_spectrometer.Manufacturer': [u'Bruker',  
u'Bruker',  
u'Bruker',  
u'Varian',  
u'Varian',  
u'Varian'],  
'_NMR_spectrometer.Model': [u'Avance',  
u'Avance',  
u'Avance',  
u'INOVA',  
u'INOVA',
```

```
u' INOVA' ]}
```

To view all of the tags in the NMR-STAR schema and their meanings, please go [here](#).

Assigned Chemical Shifts

“I just want to get the chemical shift data as an array - how do I do that?”

Keep in mind that an entry may have multiple sets of assigned chemical shifts. (For examples, there made be two sets of assignments that were made under two different sample conditions.) So to get the chemical shifts it is best to iterate through all the assigned chemical shift loops:

```
cs_result_sets = []
for chemical_shift_loop in entry15000.get_loops_by_category("
Atom_chem_shift"):
    cs_result_sets.append(chemical_shift_loop.get_tag(['Comp_
index_ID', 'Comp_ID', 'Atom_ID', 'Atom_type', 'Val', 'Val_err
']))
cs_result_sets
[[[u'2', u'SER', u'H', u'H', u'9.3070', u'0.01'],
 [u'2', u'SER', u'HA', u'H', u'4.5970', u'0.01'],
 [u'2', u'SER', u'HB2', u'H', u'4.3010', u'0.01'],
 [u'2', u'SER', u'HB3', u'H', u'4.0550', u'0.01'],
 [u'2', u'SER', u'CB', u'C', u'64.6000', u'0.1'],
 [u'2', u'SER', u'N', u'N', u'121.5800', u'0.1'],
```

```
[u'3', u'ASP', u'H', u'H', u'8.0740', u'0.01'],  
[u'3', u'ASP', u'HA', u'H', u'4.5580', u'0.01'],  
[u'3', u'ASP', u'HB2', u'H', u'2.835', u'0.01'],  
...
```

Note that we used the `get_tag()` method of the loop to only pull out the tags we were concerned with. `get_tag()` accepts an array of tags in addition to a single tag.

“But I want to access the chemical shifts as numbers, not strings!”

That is easy to do. When you first load an entry it is by default loaded with all values as strings. To instead load it such that the values match the schema, simply turn on `CONVERT_DATATYPES` prior to loading it.

```
pymrstar.CONVERT_DATATYPES = True  
ent15000 = pymrstar.Entry.from_database(15000)  
cs_result_sets = []  
for chemical_shift_loop in entry15000.get_loops_by_category("Atom_chem_shift"):  
    cs_result_sets.append(chemical_shift_loop.get_tag(['Comp_index_ID', 'Comp_ID', 'Atom_ID', 'Atom_type', 'Val', 'Val_err']))  
cs_result_sets  
[[[2, u'SER', u'H', u'H', Decimal('9.3070'), Decimal('0.01')]  
,  
 [2, u'SER', u'HA', u'H', Decimal('4.5970'), Decimal('0.01')]
```

```
],
  [2, u'SER', u'HB2', u'H', Decimal('4.3010'), Decimal('0.01')],
],
  [2, u'SER', u'HB3', u'H', Decimal('4.0550'), Decimal('0.01')],
],
  [2, u'SER', u'CB', u'C', Decimal('64.6000'), Decimal('0.1')],
],
  [2, u'SER', u'N', u'N', Decimal('121.5800'), Decimal('0.1')],
],
  [3, u'ASP', u'H', u'H', Decimal('8.0740'), Decimal('0.01')],
,
  [3, u'ASP', u'HA', u'H', Decimal('4.5580'), Decimal('0.01')],
],
  [3, u'ASP', u'HB2', u'H', Decimal('2.835'), Decimal('0.01')],
],
  [3, u'ASP', u'HB3', u'H', Decimal('2.754'), Decimal('0.01')],
],
  [3, u'ASP', u'CA', u'C', Decimal('57.6400'), Decimal('0.1')],
],
  [3, u'ASP', u'N', u'N', Decimal('121.1040'), Decimal('0.1')],
],
  ...
```

This is a great opportunity to point out that if all you want is the chemical shifts, or one or two tags, you may find it significantly easier to use the [BMRB API \(chemical shift example\)](#) to fetch that data directly and on-demand rather than dealing directly with NMR-STAR at all.

For more examples of PyNMRSTAR library usage, please look [here](#).